




MRPS5 Antibody

Product Code	CSB-PA808897
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P82675
Immunogen	Synthesized peptide derived from internal of Human MRPS5.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MRPS5 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Mitochondrial 28S ribosomal protein S5; S5mt; MRP-S5;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MRPS5

Image



Immunohistochemistry analysis of paraffin-embedded human testis tissue using MRPS5 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.