




SLC25A11 Antibody

Product Code	CSB-PA269768
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02978
Immunogen	Synthesized peptide derived from internal of Human SLC25A11.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total SLC25A11 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Mitochondrial 2-oxoglutarate/malate carrier protein; OGCP; Solute carrier family 25 member 11;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC25A11

Image



Western blot analysis of extracts from HeLa cells, using SLC25A11 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.