




MRPL3 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA553771 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P09001 |
| Immunogen | Synthesized peptide derived from Internal of Human MRPL3. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Specificity | The antibody detects endogenous levels of total MRPL3 protein. |
| Tested Applications | ELISA,WB;WB:1:500-1:3000 |
| Form | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Alias | 39S ribosomal protein L3; mitochondrial; L3mt; mitochondrial 60S ribosomal protein L3; mitochondrial ribosomal protein L3 |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | MRPL3 |

Image



Western blot analysis of extracts from K562 cells, using MRPL3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.