




SLC16A2 Antibody

Product Code	CSB-PA062013
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36021
Immunogen	Synthesized peptide derived from internal of Human SLC16A2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total SLC16A2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Monocarboxylate transporter 8; MCT 8; MCT 7; Solute carrier family 16 member 2; X-linked PEST-containing transporter
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC16A2

Image



Western blot analysis of extracts from K562 cells, using SLC16A2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.