




MEGF9 Antibody

Product Code	CSB-PA071033
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1U4
Immunogen	Synthesized peptide derived from internal of Human MEGF9.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total MEGF9 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	egf15; megf9; multiple egf-like-domains 9;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MEGF9

Image



Western blot analysis of extracts from COLO cells, using MEGF9 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.