



MYOM1 Antibody

Product Code	CSB-PA964006
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52179
Immunogen	Synthesized peptide derived from internal of Human MYOM1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MYOM1 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Myomesin-1; Myomesin family member 1; 190 kDa titin-associated protein; 190 kDa connectin-associated protein;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MYOM1

Image



Immunohistochemistry analysis of paraffin-embedded human skeletal muscle tissue using MYOM1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.