




MYH4 Antibody

Product Code	CSB-PA729938
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y623
Immunogen	Synthesized peptide derived from C-terminal of Human MYH4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MYH4 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Myosin-4; Myosin heavy chain 4; Myosin heavy chain 2b; MyHC-2b; Myosin heavy chain
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MYH4

Image



Immunohistochemistry analysis of paraffin-embedded human skeletal muscle tissue using MYH4 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.