




MYO5A Antibody

Product Code	CSB-PA785452
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4I1
Immunogen	Synthesized peptide derived from C-terminal of Human MYO5A.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MYO5A protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Myosin-Va; Dilute myosin heavy chain; non-muscle; Myosin-12; Myosin heavy chain 12
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MYO5A

Image



Western blot analysis of extracts from HeLa cells, using MYO5A antibody.