




MYO5A Antibody

Product Code	CSB-PA785452
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4I1
Immunogen	Synthesized peptide derived from C-terminal of Human MYO5A.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MYO5A protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Myosin-Va; Dilute myosin heavy chain; non-muscle; Myosin-12; Myosin heavy chain 12
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MYO5A

Image



Western blot analysis of extracts from HeLa cells, using MYO5A antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.