




NDUFA4 Antibody

Product Code	CSB-PA955890
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00483
Immunogen	Synthesized peptide derived from internal of Human NDUFA4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total NDUFA4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CI-MLRQ; complex I-MLRQ; EC 1.6.5.3; EC 1.6.99.3; MLRQ
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFA4

Image



Western blot analysis of extracts from HepG2 cells and RAW264.7 cells, using NDUFA4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.