




NDUFA4L2 Antibody

Product Code	CSB-PA563041
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRX3
Immunogen	Synthesized peptide derived from internal of Human NDUFA4L2.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total NDUFA4L2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 4-like 2; NADH-ubiquinone oxidoreductase MLRQ subunit homolog; NUOMS; NDUFA4L2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFA4L2

Image



Western blot analysis of extracts from Jurkat cells, COLO cells, HUVEC cells and HepG2 cells, using NDUFA4L2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.