




NDUFB1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA567533 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O75438 |
| Immunogen | Synthesized peptide derived from internal of Human NDUFB1. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Specificity | The antibody detects endogenous levels of total NDUFB1 protein. |
| Tested Applications | ELISA,IHC;IHC:1:50-1:100 |
| Form | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Alias | NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 1; NADH-ubiquinone oxidoreductase MNLL subunit; Complex I-MNLL; CI-MNLL; |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | NDUFB1 |

Image



Immunohistochemistry analysis of paraffin-embedded human pancreas tissue using NDUFB1 antibody.