




NDUFB9 Antibody

Product Code	CSB-PA783045
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6M9
Immunogen	Synthesized peptide derived from internal of Human NDUFB9.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total NDUFB9 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	B22; CI-B22; complex I-B22; EC 1.6.5.3; EC 1.6.99.3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFB9

Image



Western blot analysis of extracts from COLO cells and 293 cells, using NDUFB9 antibody.