




NDUFV3 Antibody

Product Code	CSB-PA049345
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56181
Immunogen	Synthesized peptide derived from internal of Human NDUFV3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total NDUFV3 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	NADH dehydrogenase [ubiquinone] flavoprotein 3; mitochondrial; NADH-ubiquinone oxidoreductase 9 kDa subunit; Complex I-9kD; CI-9kD
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFV3

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue using NDUFV3 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.