





NDUFS5 Antibody

Product Code	CSB-PA942514
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43920
Immunogen	Synthesized peptide derived from internal of Human NDUFS5.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total NDUFS5 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	NADH dehydrogenase [ubiquinone] iron-sulfur protein 5; NADH-ubiquinone oxidoreductase 15 kDa subunit; Complex I-15 kDa; CI-15 kDa;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFS5

Image



Western blot analysis of extracts from K562 cells, using NDUFS5 antibody.



Immunohistochemistry analysis of paraffin-embedded human heart tissue using NDUFS5 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.