




CYB5R1 Antibody

Product Code	CSB-PA445001
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHQ9
Immunogen	Synthesized peptide derived from internal of Human CYB5R1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total CYB5R1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	b5R.1; cytochrome b5 reductase 1; cytochrome b5 reductase 1 (B5R.1); EC 1.6.2.2; H:quinone oxidoreductase type 3 polypeptide A2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYB5R1

Image



Western blot analysis of extracts from RAW264.7 cells, using CYB5R1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.