





CYB5R3 Antibody

Product Code	CSB-PA178677
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00387
Immunogen	Synthesized peptide derived from internal of Human CYB5R3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CYB5R3 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	B5R; cytochrome b5 reductase; cytochrome b5 reductase 3; DIA1; diaphorase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYB5R3

Image



Western blot analysis of extracts from COLO cells, HUVEC cells, HT-29 cells, K562 cells and HepG2 cells, using CYB5R3 antibody.



Immunohistochemistry analysis of paraffin-embedded human placenta tissue using CYB5R3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.