





MT-ND1 Antibody

Product Code	CSB-PA267250
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03886
Immunogen	Synthesized peptide derived from internal of Human MT-ND1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MT-ND1 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	NADH-ubiquinone oxidoreductase chain 1; EC 1.6.5.3; NADH dehydrogenase subunit 1; MT-ND1; MTND1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MT-ND1

Image



Western blot analysis of extracts from A549 cells and COLO cells, using MT-ND1 antibody.



Immunohistochemistry analysis of paraffin-embedded human placenta tissue, using MT-ND1 antibody.