





# MT-ND1 Antibody

<b>Product Code</b>	CSB-PA267250
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P03886
<b>Immunogen</b>	Synthesized peptide derived from internal of Human MT-ND1.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Specificity</b>	The antibody detects endogenous levels of total MT-ND1 protein.
<b>Tested Applications</b>	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
<b>Form</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	NADH-ubiquinone oxidoreductase chain 1; EC 1.6.5.3; NADH dehydrogenase subunit 1; MT-ND1; MTND1
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	MT-ND1

## Image



Western blot analysis of extracts from A549 cells and COLO cells, using MT-ND1 antibody.



Immunohistochemistry analysis of paraffin-embedded human placenta tissue, using MT-ND1 antibody.

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.