




MT-ND5 Antibody

Product Code	CSB-PA680236
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03915
Immunogen	Synthesized peptide derived from internal of Human MT-ND5.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MT-ND5 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	NADH-ubiquinone oxidoreductase chain 5; EC=1.6.5.3; NADH dehydrogenase subunit 5; MTND5; NADH5
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MT-ND5

Image



Western blot analysis of extracts from HT-29 cells, using MT-ND5 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.