




NOX1 Antibody

Product Code	CSB-PA924751
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5S8
Immunogen	Synthesized peptide derived from C-terminal of Human NOX1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total NOX1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	NADPH oxidase 1; NOX-1; NOH-1; NADH/NADPH mitogenic oxidase subunit P65-MOX; Mitogenic oxidase 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NOX1

Image



Western blot analysis of extracts from 293 cells, using NOX1 antibody.