




NOX1 Antibody

Product Code	CSB-PA924751
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5S8
Immunogen	Synthesized peptide derived from C-terminal of Human NOX1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total NOX1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	NADPH oxidase 1; NOX-1; NOH-1; NADH/NADPH mitogenic oxidase subunit P65-MOX; Mitogenic oxidase 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NOX1

Image



Western blot analysis of extracts from 293 cells, using NOX1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.