




N4BP1 Antibody

Product Code	CSB-PA620222
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75113
Immunogen	Synthesized peptide derived from internal of Human N4BP1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total N4BP1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Nedd4 binding protein 1; NEDD4-binding protein 1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	N4BP1

Image



Western blot analysis of extracts from MCF-7 cells, 293 cells and Jurkat cells, using N4BP1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.