




NPHP4 Antibody

Product Code	CSB-PA870638
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75161
Immunogen	Synthesized peptide derived from internal of Human NPHP4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total NPHP4 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Nephrocystin-4; Nephroretinin;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NPHP4

Image



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using NPHP4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.