




NRSN1 Antibody

Product Code	CSB-PA967502
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZ57
Immunogen	Synthesized peptide derived from internal of Human NRSN1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total NRSN1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	neurensin 1; nrsn1; p24; vesicular membrane protein p24; vmp
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NRSN1

Image



Western blot analysis of extracts from HeLa cells, using NRSN1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.