





NXPH4 Antibody

Product Code	CSB-PA448327
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95158
Immunogen	Synthesized peptide derived from internal of Human NXPH4.
Raised In	Rabbit
Species Reactivity	Human,Rat
Specificity	The antibody detects endogenous levels of total NXPH4 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	neurexophilin-4; NPH4; NXPH4;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NXPH4

Image



Western blot analysis of extracts from COLO cells and HUVEC cells, using NXPH4 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using NXP4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.