




NBPF1/NBPF9/NBPF10/NBPF12/NBPF14/NBPF15/NBPF16/NBPF20 Antibody

Product Code	CSB-PA279287
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TAG4
Immunogen	Synthesized peptide derived from internal of Human NBPF1/9/10/12/14/15/16/20.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total NBPF1/9/10/12/14/15/16/20 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	NBPF1; neuroblastoma breakpoint family; 12;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NBPF12/NBPF10/NBPF1/NBPF9

Image



Western blot analysis of extracts from JK cells, using NBPF1/9/10/12/14/15/16/20 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.