




NBPF7 Antibody

Product Code	CSB-PA052635
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C2Y1
Immunogen	Synthesized peptide derived from internal of Human NBPF7.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total NBPF7 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Neuroblastoma breakpoint family member 7; NBPF7;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NBPF7

Image



Western blot analysis of extracts from HepG2 cells and HT-29 cells, using NBPF7 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.