




NPAS4 Antibody

Product Code	CSB-PA159106
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUM7
Immunogen	Synthesized peptide derived from internal of Human NPAS4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total NPAS4 protein.
Tested Applications	ELISA, IHC; IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Neuronal PAS domain-containing protein 4; Neuronal PAS4; HLH-PAS transcription factor NXF;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NPAS4

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue using NPAS4 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.