





VGF Antibody

Product Code	CSB-PA220007
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15240
Immunogen	Synthesized peptide derived from internal of Human VGF.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total VGF protein
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	neuro-endocrine specific VGF; neurosecretory protein VGF; VGF nerve growth factor inducible;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	VGF

Image



Western blot analysis of extracts from RAW264.7 cells, using VGF antibody.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using VEGF antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.