




NCOA7 Antibody

Product Code	CSB-PA127723
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NI08
Immunogen	Synthesized peptide derived from internal of Human NCOA7.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total NCOA7 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ERAP140; ESNA1; estrogen nuclear receptor coactivator 1; estrogen receptor associated protein 140 kDa; nuclear receptor coactivator 7
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NCOA7

Image



Western blot analysis of extracts from LOVO cells, 3T3 cells and RAW264.7 cells, using NCOA7 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.