




ZDHHC15 Antibody

Product Code	CSB-PA295024
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MV8
Immunogen	Synthesized peptide derived from C-terminal of Human ZDHHC15.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Specificity	The antibody detects endogenous levels of total ZDHHC15 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Palmitoyltransferase ZDHHC15; EC 2.3.1.-; Zinc finger DHHC domain-containing protein 15; DHHC-15; ZDHHC15
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZDHHC15

Image



Western blot analysis of extracts from Jurkat cells, using ZDHHC15 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.