




ZDHHC7 Antibody

Product Code	CSB-PA992709
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NXF8
Immunogen	Synthesized peptide derived from C-terminal of Human ZDHHC7.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZDHHC7

Image



Western blot analysis of extracts from HUVEC cells, using ZDHHC7 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.