





EHHADH Antibody

Product Code	CSB-PA094370
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08426
Immunogen	Synthesized peptide derived from internal of Human EHHADH.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total EHHADH protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 1.1.1.35; EC 4.2.1.17; EC 5.3.3.8; ECHD; ECHP
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EHHADH

Image



Western blot analysis of extracts from A549 cells, using EHHADH antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using EHHADH antibody.