





PEX14 Antibody

Product Code	CSB-PA775947
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75381
Immunogen	Synthesized peptide derived from internal of Human PEX14.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total PEX14 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Peroxin-14; Peroxisomal membrane anchor protein PEX14; Peroxisomal membrane protein PEX14; PEXE; PTS1 receptor docking protein
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PEX14

Image



Western blot analysis of extracts from A549 cells, using PEX14 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using PEX14 antibody.