





PEX7 Antibody

Product Code	CSB-PA114293
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00628
Immunogen	Synthesized peptide derived from internal of Human PEX7.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PEX7 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Peroxisomal targeting signal 2 receptor; PTS2 receptor; Peroxin-7;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PEX7

Image



Western blot analysis of extracts from NIH/3T3 cells, using PEX7 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using PEX7 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.