




PEBP1 Antibody

Product Code	CSB-PA948956
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30086
Immunogen	Synthesized peptide derived from internal of Human PEBP1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PEBP1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	23 kDa morphine-binding protein; B-RAF proto-oncogene serine/threonine-protein kinase; HCNP; Hippocampal cholinergic neurostimulating peptide; Neuropolypeptide h3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PEBP1

Image



Western blot analysis of extracts from K562 cells and Jurkat cells, using PEBP1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.