




PITPNB Antibody

Product Code	CSB-PA877125
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48739
Immunogen	Synthesized peptide derived from internal of Human PITPNB.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total PITPNB protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Phosphatidylinositol transfer protein beta isoform; PI-TP-beta; PI-TPbeta; PIPNB; PPI2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PITPNB

Image



Western blot analysis of extracts from HepG2 cells, using PITPNB antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.