




PIK3R5 Antibody

Product Code	CSB-PA288007
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WYR1
Immunogen	Synthesized peptide derived from internal of Human PIK3R5.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total PIK3R5 protein.
Tested Applications	ELISA, WB, IHC, IF; WB: 1:500-1:3000, IHC: 1:50-1:200, IF: 1:100-1:200
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Phosphoinositide 3-kinase regulatory subunit 5; PI3-kinase regulatory subunit 5; PI3-kinase p101 subunit; PtdIns-3-kinase p101; p101-PI3K
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PIK3R5

Image



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using PIK3R5 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.