




PLD4 Antibody

Product Code	CSB-PA219574
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BZ4
Immunogen	Synthesized peptide derived from internal of Human PLD4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total PLD4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	C14orf175; choline phosphatase 4; EC 3.1.4.4; phosphatidylcholine-hydrolyzing phospholipase D4; phospholipase D4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PLD4

Image



Western blot analysis of extracts from Jurkat cells, using PLD4 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.