




KCNA5 Antibody

Product Code	CSB-PA074982
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22460
Immunogen	Synthesized peptide derived from internal of Human KCNA5.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total KCNA5 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Potassium voltage-gated channel subfamily A member 5; Voltage-gated potassium channel subunit Kv1.5; HK2; HPCN1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KCNA5

Image



Western blot analysis of extracts from 293 cells, using KCNA5 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.