




# KCNQ4 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA288698   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | P56696   |
| <b>Immunogen</b>           | Synthesized peptide derived from C-terminal of Human KCNQ4.  |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Specificity</b>         | The antibody detects endogenous levels of total KCNQ4 protein.   |
| <b>Tested Applications</b> | ELISA, WB; WB: 1:500-1:3000  |
| <b>Form</b>                | Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.         |
| <b>Purification Method</b> | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                      |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Potassium voltage-gated channel subfamily KQT member 4; Voltage-gated potassium channel subunit Kv7.4; Potassium channel subunit alpha KvLQT4; KQT-like 4; |
| <b>Product Type</b>        | Polyclonal Antibody  |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Target Names</b>        | KCNQ4  |

## Image



Western blot analysis of extracts from HuvEc cells, using KCNQ4 antibody.

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.