





NDRG4 Antibody

Product Code	CSB-PA189877
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULP0
Immunogen	Synthesized peptide derived from internal of Human NDRG4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total NDRG4 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Brain development-related molecule 1; Protein NDRG4; SMAP-8; Vascular smooth muscle cell-associated protein 8;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NDRG4

Image



Western blot analysis of extracts from HepG2 cells and HT-29 cells, using NDRG4 antibody.



Immunohistochemistry analysis of paraffin-embedded human heart tissue using NDRG4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.