




PPP1R7 Antibody

Product Code	CSB-PA788208
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15435
Immunogen	Synthesized peptide derived from internal of Human PPP1R7.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PPP1R7 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	PP1R7; Protein phosphatase 1; regulatory (Inhibitor) subunit 7; Protein phosphatase 1; regulatory subunit 7
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPP1R7

Image



Western blot analysis of extracts from HUVEC cells, using PPP1R7 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.