




S100A16 Antibody

Product Code	CSB-PA113412
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96FQ6
Immunogen	Synthesized peptide derived from internal of Human S100A16.
Raised In	Rabbit
Species Reactivity	Human,Rat
Specificity	The antibody detects endogenous levels of total S100A16 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	S100 calcium-binding protein A16; S100F; S10AG; S116;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	S100A16

Image



Western blot analysis of extracts from 3T3 cells and RAW264.7 cells, using S100A16 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.