




S100Z Antibody

Product Code	CSB-PA617327
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXG8
Immunogen	Synthesized peptide derived from internal of Human S100Z.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total S100Z protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Protein S100-Z; S100 calcium-binding protein Z;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	S100Z

Image



Western blot analysis of extracts from COLO cells, using S100Z antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.