





ABCC13 Antibody

Product Code	CSB-PA870553
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NSE7
Immunogen	Synthesized peptide derived from internal of Human ABCC13.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ABCC13 protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ABCCD; ATP-binding cassette protein C13; PRED6;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABCC13

Image



Western blot analysis of extracts from 293 cells, COLO cells and HT-29 cells, using ABCC13 antibody.



Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue using ABCC13 antibody.