




APOBEC4 Antibody

Product Code	CSB-PA578412
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WW27
Immunogen	Synthesized peptide derived from internal of Human APOBEC4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total APOBEC4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	utative C->U-editing enzyme APOBEC-4; EC 3.5.4.-; Apolipoprotein B mRNA-editing enzyme catalytic polypeptide-like 4; APOBEC4; C1orf169
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	APOBEC4

Image



Western blot analysis of extracts from K562 cells, using APOBEC4 antibody.