




PWWP2B Antibody

Product Code	CSB-PA149091
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NUJ5
Immunogen	Synthesized peptide derived from internal of Human PWWP2B.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PWWP2B protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	PWP2B; PWWP domain-containing protein 2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PWWP2B

Image



Western blot analysis of extracts from LOVO cells, using PWWP2B antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.