




PDK2 Antibody

Product Code	CSB-PA574849
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15119
Immunogen	Synthesized peptide derived from internal of Human PDK2.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total PDK2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.11.2; kinase Pyruvate dehydrogenase kinase 2; PDK2; Pyruvate dehydrogenase [lipoamide] kinase isozyme 2; mitochondrial precursor
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PDK2

Image



Western blot analysis of extracts from HepG2 cells, using PDK2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.